

medilogTMFD5

DIGITAL HOLTER RECORDER

Every beat of the heart... every second of the day

The MedilogFD5 is the next step in the evolution of digital Holter recorders. With a simple user interface and compact design, it provides all the capabilities required for reliable ECG recording.

- View ECG at hook-up
- Simple operation
- Records for up to 48 hours
- Compact design



Specification

Input channels:	3 channel recording mode
Recording duration:	24 or 48 hours
Dimensions:	70 x 110 x 20mm (2.75" x 4.37" x 0.80")
Weight with batteries:	141g (5oz)
Display:	LCD for ECG and recorder setup
Storing rate:	128 Hz for each channel
Sampling resolution:	10 bit
Gain:	Adjustable 0.5x, 1x and 2x
Temperature:	0 to 45°C
Humidity:	10 to 95% non-condensing
Power:	2 x AA type alkaline batteries

View ECG at hook-up

The liquid crystal display enables verification of high quality ECG signal at hook-up. This allows remote patient connection and eliminates the need for separate signal monitoring hardware. In addition, viewing high quality ECG at hook-up ensures a better quality ECG and enables more efficient analysis.

Simple operation

Simple keypad and on-screen controls make changing settings and recording an easy task. All functions can be accessed quickly and easily allowing you to spend more time with your patients.

Standard USB connection allows easy download of ECG to the Medilog Holter Management System without the need for additional hardware.

Records up to 48 hours

The Medilog FD5 can record for up to 48 hours. This increases the opportunity to capture significant clinical events and meets the need for extended recordings.

Compact design

Compact and lightweight, the Medilog FD5 combines the rigours of ambulatory monitoring with the need for patient comfort. This inherently improves patient compliance for long-term ambulatory monitoring.

The sleek, practical designed recorder has a removable input lead. The same connection is used for the fast download of ECG recordings.

This product complies with the EN 60601-1 safety standard. The CE mark indicates that it meets the essential requirements of European Council Directive 93/42/EEC of 14 June 1993 concerning medical devices. Huntleigh Healthcare operates a quality system certified to ISO 9001:2000, ISO 13485:2003, and the TickIT Guide Issue 5. Medilog is a registered trademark of Huntleigh Healthcare. Microsoft and Windows are registered trademarks of Microsoft Corporation. For the USA the following applies: CAUTION: Federal law restricts this device to sale by or on the order of a physician.

Copyright © 2005 Huntleigh Healthcare. All rights reserved.
Huntleigh Healthcare, Diagnostic Products Division, 35 Portmanmoor Road, Cardiff CF24 5HN, UK.



0088



Part Number: APL535/06/05/EN

CARDIOLOGY Products Division

Lion House, Oriental Road
Woking, Surrey
GU22 8AP, United Kingdom

CARDIOLOGY Products Division

Otto-von-Guericke-Ring 10
65205 Wiesbaden, Germany
T: +49 (0) 6122 937 150
F: +49 (0) 6122 937 148

www.medilogholter.com

 **Huntleigh**
HEALTHCARE